

# cardioWell

Cardiowell's 'always-connected' cellular-enabled blood pressure monitor automatically sends results thru the secure Cardiowell Cloud Network, allowing results to be immediately available on the Cardiowell App.

The Cardiowell 'always-connected' blood pressure monitor works straight out of the box.

Easy-to-use, one-button operation provides accurate readings in under 30 seconds. Captures blood pressure readings on inflation, without over inflating or causing discomfort.

Interested in using the device for a study today? Get on the waiting list for the FDA cleared release.

Please contact [info@welladapt.com](mailto:info@welladapt.com) for details.

## ALWAYS CONNECTED BLOOD PRESSURE MONITOR



### User-friendly

Simple setup. Open the box, press start and view results



### Clinical Accuracy

Clinical grade results in under 30 seconds.

### Health Management

Taking and tracking blood pressure can improve your health.



### Connectivity

Cellular technology enables real-time blood pressure tracking anywhere you go.

Investigational Use Only. Pending FDA Approval

Learn more at  
[cardiowell.io](http://cardiowell.io)

315 Meigs Rd Suite A385, Santa Barbara, CA 93109, USA

[www.cardiowell.io](http://www.cardiowell.io)

[info@welladapt.com](mailto:info@welladapt.com)

# cardioWell

B L O O D P R E S S U R E M O N I T O R

## Product Details

### User-friendly Design:

Big button and LCD display for users to operate easier.  
Measurement of Pressure and Pulse

### SID (Simultaneously Information Design):

User can read systolic, diastolic, and pulse at the same time.

### Power-Saving Design:

Auto power off after two (2) minutes.

### Low Battery indication:

Low battery symbol is shown when battery low power

### Memory recall:

Users can recall up to 60 measurements with date/ time stamp, plus an average of the last 3 stored measurements in memory to share with one's physician or a trained healthcare professional.

**Alarm Clock (Optional)** The alarm clock can be set to ring at three different times in the day. The alarm clock is especially helpful if users are required to measure blood pressure at specific times.

### IHB Indicator (Optional)

Automatically detect IHB and show the IHB icon if IHB is detected when finishing measurement.

### WHO Indicator (Optional):

The WHO Indicator, located on the bottom of LCD Display. It will indicate your reading in a range (normal, high-normal, high) in accordance with standards set by the World Health Organization (WHO).

**Device Monitoring:** Transmits blood pressure, battery level and cellular signal strength over the AT&T GSM network to the CardioWell cloud network allowing remote troubleshooting of the device.



CardioWell

Download from the  
Apple App Store

## Technical Specifications

**Measuring Range:** Blood Pressure: 30 – 280 mmHg;

**Pulse Rate:** 40 – 199 beats/min

**Calibration Accuracy:** Blood Pressure: +/- 3 mmHg; Pulse: +/- 4% of reading

**Operating Environment:** -10°C to 40°C (50°F to 104°F), with relative humidity up to 85% (non-condensing)

**Storage Environment:** -20°C to 50°C (-4°F to 122°F), with relative humidity up to 85% (non-condensing)

**Power Source:** 4 x 1.5 V AA alkaline (included)

**Weight:** 492 g / 1.1 lb, excluding batteries

**Dimensions:** 130 mm x 150 mm x 65 mm (W x H x D) ; 5.1" x 5.9" x 2.6" (W x H x D)

**Battery Life:** 365 measurements

**Cuff Circumference:** 229 mm – 330 mm (included); 9" – 13"

**Connectivity:** GSM

**Accessories:** AC Adaptor, S, M, L cuffs

*Specifications subject to change*



Configurable blood pressure reminders

Learn more at  
**cardiowell.io**

315 Meigs Rd Suite A385, Santa Barbara, CA 93109, USA

 [www.cardiowell.io](http://www.cardiowell.io)

 [info@welladapt.com](mailto:info@welladapt.com)